


## HighFidelity 16 PCR Labeling Kit (Biotin)

Cat #: orb653771 (manual)



Structural formula of HighFidelity Biotin16 PCR Labeling Kit

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

**Additional Storage Conditions:** avoid freeze/thaw cycles

**Shelf Life:** 12 months

### Description:

HighFidelity Biotin16 PCR Labeling Kit is designed to produce randomly Biotin16-modified DNA probes by PCR. Such probes are ideally suited for *in situ* hybridization and Northern Blot experiments. PCR-based labeling is superior to random-primed labeling by Klenow fragment if template amounts are limited or amplification of a specific DNA fragments is required.

Biotin-16-dUTP is efficiently incorporated into DNA as substitute for its natural counterpart dTTP using an optimized reaction buffer and a High Fidelity Polymerase. 50 % Biotin-16-dUTP substitution typically results in an optimal balance between reaction and labeling efficiency. Individual optimization of Biotin-16-dUTP/dTTP ratio however, can easily be achieved with the single nucleotide format. The resulting Biotin16-modified DNA probe can subsequently be detected by HRP- or AP-modified Streptavidin.

The kit contains sufficient reagents for 175 labeling reactions of 20 µl each (50 % Biotin16-dUTP substitution, 100 µM dATP/dGTP/dCTP, 50 µM dTTP, 50 µM Biotin16-dUTP).

### Content:

#### High Fidelity Polymerase

in storage buffer with 50% glycerol (v/v) 1x 200  $\mu$ l (500 units, 2.5 units/ $\mu$ l)

### **High Fidelity Labeling Buffer**

1x 500  $\mu$ l (10x)

#### **dATP - Solution**

1x 20  $\mu$ l (100 mM)

#### **dGTP - Solution**

1x 20  $\mu$ l (100 mM)

#### **dCTP - Solution**

1x 20  $\mu$ l (100 mM)

#### **dTTP - Solution**

1x 20  $\mu$ l (100 mM)

#### **Biotin-16-dUTP**

1x 200  $\mu$ l (1 mM)

#### **Lambda DNA**

1x 20  $\mu$ l (100 ng/ $\mu$ l)

#### **500 bp forward primer**

1x 20  $\mu$ l (10  $\mu$ M)

#### **500 bp reverse primer**

1x 20  $\mu$ l (10  $\mu$ M)

#### **PCR-grade water**

1x 1.2 ml

#### **To be provided by user**

DNA template Primer

DNA purification tools (optional)

## **1. Preparation of working solutions**

### 1.1 Preparation of 1 mM dATP/dCTP/dGTP working solution

- Thaw 100 mM dATP, 100 mM dCTP and 100 mM dGTP solutions on ice, vortex and spin-down briefly.
- Prepare a 1:100 dilution with PCR-grade water to achieve a final concentration of 1 mM (e.g. 2  $\mu$ l 100 mM dATP + 2  $\mu$ l 100 mM dCTP + 2  $\mu$ l 100 mM dGTP + 194  $\mu$ l PCR-grade water).
- 1 mM ATP/CTP/GTP working solution can be stored at  $-20^{\circ}\text{C}$ . Pre- pare aliquots to avoid freeze/thaw cycles.

### 1.2 Preparation of 1 mM dTTP working solution

- Thaw 100 mM dTTP solution on ice, vortex and spin-down briefly.

Explore. Bioreagents.

- Prepare a 1:100 dilution with PCR-grade water to achieve a final concentration of 1 mM (e.g. 2  $\mu$ l 100 mM dTTP + 198  $\mu$ l PCR-grade water).
- 1 mM dTTP working solution can be stored at -20 °C. Prepare aliquots to avoid freeze/thaw cycles.

## 2. Standard PCR Labeling protocol

The standard protocol is set-up for labeling of a 500 bp DNA fragment. An optimal balance between reaction and labeling efficiency is typically achieved with 50% Biotin-16-dUTP substitution following the standard protocol below however, individual optimization might improve results for individual applications.

- Assemble the PCR on ice in the order stated below (DNase-free reaction tube).
- Vortex and spin-down briefly.

Component	Volume	Final concentration
PCR-grade water	X $\mu$ l	
High Fidelity Buffer (10x)	2 $\mu$ l	1x
1 mM dATP/dCTP/ dGTP working solution (s. 1.1)	2 $\mu$ l	100 $\mu$ M
1 mM dTTP working solution (s. 1.2)	1 $\mu$ l	50 $\mu$ M
1 mM Biotin-16- dUTP	1 $\mu$ l	50 $\mu$ M
forward primer (10 $\mu$ M)	X $\mu$ l	0.1 - 1 $\mu$ M (e.g. 0.3 $\mu$ M 500 bp forward primer)
reverse primer (10 $\mu$ M)	X $\mu$ l	0.1 - 1 $\mu$ M (e.g. 0.3 $\mu$ M 500 bp reverse primer)
template DNA	X $\mu$ l	1 - 10 ng genomic DNA (e.g. 1 ng Lambda DNA)
High Fidelity Pol (2.5 units/ $\mu$ l)	1 $\mu$ l	2.5 units
Total volume	20 $\mu$ l	

## Recommended cycling conditions

Cycle step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1x
Denaturation	95°C	20 sec	30x
Annealing1)	58°C	30 sec	
Elongation2)	68°C	60 sec	
Final Elongation	68°C	2 min	1x

- 1)The annealing temperature depends on the melting temperature of primers used.
- 2)The elongation time depends on the length of fragments to be amplified. A time of 2 min/kbp is recommended. Elongation at 72°C works as well.

For optimal amplification results and high incorporation rates an individual optimization of the recommended PCR assay and cycling conditions may be necessary for each new primer-template pair.

### **3. Probe purification:**

Probe purification is not required for most hybridization experiments. If a downstream application requires purification (e.g. concentration determination by absorbance measurement) we recommend silica-membrane or gel filtration-based purification.