

**HoxD10 rabbit pAb****Cat#: orb768642 (Manual)**


For research use only. Not intended for diagnostic use.

<b>Product Name</b>	HoxD10 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB;ELISA
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human HOXD10. AA range:291-340
<b>Specificity</b>	HoxD10 Polyclonal Antibody detects endogenous levels of HoxD10 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Homeobox protein Hox-D10
<b>Gene Name</b>	HOXD10
<b>Cellular localization</b>	Nucleus.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal

<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	32kD
<b>Human Gene ID</b>	3236
<b>Human Swiss-Prot Number</b>	P28358
<b>Alternative Names</b>	HOXD10; HOX4D; HOX4E; Homeobox protein Hox-D10; Homeobox protein Hox-4D; Homeobox protein Hox-4E

**Background**

This gene is a member of the Abd-B homeobox family and encodes a protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox D genes located on chromosome 2. The encoded nuclear protein functions as a sequence-specific transcription factor that is expressed in the developing limb buds and is involved in differentiation and limb development. Mutations in this gene have been associated with Wilm's tumor and congenital vertical talus (also known as "rocker-bottom foot" deformity or congenital convex pes valgus) and/or a foot deformity resembling that seen in Charcot-Marie-Tooth disease. [provided by RefSeq, Jul 2008],



**Western blot analysis of lysates from K562 and RAW264.7 cells, using HOXD10 Antibody. The lane on the right is blocked with the synthesized peptide.**