

Product Datasheet

Vazyme - 5min Universal Ligation Mix (C311-02)

Description 5min Universal Ligation Mix is a ready-to-use 2 x premix solution, containing T4 DNA ligase that ca...

Storage Store at -30 ~ -15°C and transport at ≤0°C.

Note For research use only

Application notes Compatible: Compatible with sticky ends, blunt ends, TA , linker or adapter ligation Fast: Ligation can be completed quickly in 5 minutes at 25°C Convenient: Reaction mixture can be used directly for bacterial transformation

Expiration Date 12 months from date of receipt.

