

Product Datasheet

Vazyme - FastPure Blood/Cell/Tissue/Bacteria DNA Isolation Mini Kit (DC112-02)

Description

This kit is suitable for extracting genomic DNA from $\leq 200 \mu\text{l}$ of fresh or frozen anticoagulated wh...

**Storage**

Store at 15 ~ 25°C and transport at room temperature. If ambient temperatures often exceed 25°C, we suggest storing Proteinase K at 2 ~ 8°C.

Note

For research use only

Application notes

High quality and good integrity: The extracted DNA has the characteristics of high yield, good purity and complete fragments
Wide sample applicability: Suitable for direct extraction of DNA from fresh or frozen anticoagulated blood, cells, animal tissues and bacteria samples
Simple and efficient: Extraction from animal tissue can be completed rapidly in 30 min

Expiration Date

12 months from date of receipt.