

Product Datasheet

Vazyme - FastPure Gel DNA Extraction Mini Kit (DC301-01)

Description	This kit adopts optimized buffer system and silica gel column purification technology, which can rec...
Storage	Store at 15 ~ 25°C and transport at room temperature.
Note	For research use only
Application notes	One kit for two uses: Applicable to the recovery of DNA from gel and purification of crude DNA products Simple and fast: Just add 1× volume of gel solubilization buffer; no need for column equilibration; simple and fast purification in 10 - 15 min Highly efficient recovery: Up to 80% recovery efficiency Wide recovery range of DNA fragments: Compatible with DNA fragments of 70 bp - 20 kb
Expiration Date	12 months from date of receipt.

