

Product Datasheet

Vazyme - VAMNE MagUltra Blood Genomic DNA Extraction Kit (DM101-02)

Description	VAMNE MagUltra Blood Genomic DNA Extraction Kit is intended for fast and easy genomic DNA (gDNA) iso...
Storage	Store at 15 ~ 25°C and transport at room temperature. If ambient temperatures often exceed 25°C, we suggest storing Proteinase K and MagUltra Beads at 2 ~ 8°C.
Note	For research use only
Application notes	Flexible and Sensitive: Compatible with 50 µl - 1 ml EDTA, sodium citrate and heparin sodium peripheral blood, umbilical cord blood, bone marrow blood, concentrated blood and cells, etc. Fast and Convenient: Small volume extraction can be completed within 40 min by magnetic adsorption without centrifugation. High Quality: Obtain gDNA products with high yield, good purity and integrity.
Expiration Date	Please enquire.

