

Product Datasheet

Vazyme - VAHTS Target Capture Universal Blockers and Post-PCR Primer Mix for Illumina-TS (NC101-02)

Description	VAHTS Target Capture Universal Blockers and Post-PCR Primer Mix for Illumina-TS contains Universal B...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Universal adapter blocking sequence developed for common DNA libraries on the Illumina platform. Suitable for single-end or paired-end 6 - 8 bp Barcode-labeled DNA libraries.
Expiration Date	12 months from date of receipt.

