

Product Datasheet

Vazyme - 2 × Taq Plus Master Mi (P211-01)

Description	2 × Taq Plus Master Mix contains Taq Plus DNA Polymerase, dNTPs, and an optimized buffer system, wh...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Good performance: ~6X higher than Taq DNA Polymerase; higher amplification yield Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination Additional features: 3'→5' exonuclease activity Good Stability: Stable performance after 50 freeze-thaw cycles
Expiration Date	12 months from date of receipt.

