

Product Datasheet

Vazyme - 2 × Taq Plus Master Mix (P211-03)

Description 2 × Taq Plus Master Mix contains Taq Plus DNA Polymerase, dNTPs, and an optimized buffer system, wh...

Storage Store at -30 ~ -15°C and transport at ≤0°C.

Note For research use only

Application notes Good performance: ~6X higher than Taq DNA Polymerase; higher amplification yield
Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination
Additional features: 3'→5' exonuclease activity
Good Stability: Stable performance after 50 freeze-thaw cycles

Expiration Date 12 months from date of receipt.

