

## Product Datasheet

# Vazyme - 2 × Taq Plus Master Mix II (Dye Plus) (P213-01)

**Description** 2 × Taq Plus Master Mix II (Dye Plus) contains Taq Plus DNA Polymerase, dNTP and an optimized buffer...



**Storage** Store at -30 ~ -15°C for up to 12 months, and transport at ≤0°C.

**Note** For research use only

**Application notes** Add 3'→5' exonuclease activity; Fidelity is 6 times that of Taq DNA Polymerase; Improve amplification yield and length; Ready-to-use premix, minimizing experimental steps

**Expiration Date** 12 months from date of receipt.