

Product Datasheet

Vazyme - 2 × Rapid Taq Master Mix (P222-04)

Description 2 × Rapid Taq Master Mix contains Taq DNA Polymerase, elongation promoting factor, dNTP and an opti...

Storage Store at -30 ~ -15°C and transport at ≤0°C. It can be stored at 4°C for 3 months after thawing.

Note For research use only

Application notes Rapid amplification: 15 sec/kb; extension speed up to 1 sec/kb (≤ 1kb fragments)
Convenience: Include tracking dye, direct loading of PCR product on gel
Good stability: Stable performance after 50 freeze-thaw cycles
Ready-to-use: Simplified workflow and minimized pipetting steps with 2 × master mix format

Expiration Date 12 months from date of receipt.

