

## Product Datasheet

### Vazyme - 2 × Vazyme LAmp Master Mix (P311-03)

<b>Description</b>	Vazyme LAmp DNA Polymerase is a mixed enzyme of Taq DNA Polymerase and a protein containing 3'→5' ...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Good performance: ~6X higher than Taq DNA Polymerase; higher amplification efficiency Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination Additional features: 3'→5' exonuclease activity Broad Compatibility: Amplify fragments up to 21 kb (gDNA)
<b>Expiration Date</b>	12 months from date of receipt.

