

Product Datasheet

Vazyme - ChamQ SYBR qPCR Master Mix (High ROX Premixed) (Q341-02)

Description	ChamQ SYBR qPCR Master Mix (High ROX Premixed) is a special premix for qPCR assays based on SYBR Gre...
Storage	Store at -30 ~ -15°C and protect from light. Transport at ≤0°C. After thawing, the Master mix can be stored stably for 6 months at 2 ~ 8°C and protected from light.
Note	For research use only
Application notes	High efficiency hot-start polymerase: Novel antibody modified hot-start DNA polymerase with strong amplification performance and high sensitivity. High specificity: The specificity-enhancing factor Exactor inhibits non-specific amplification and primer dimerization to a large extent. High concentration ROX is premixed in Biorbyt #Q331. It is applicable for Step One, Step One Plus, ABI 5700/7300/7700/7900/7900HT and other qPCR instruments.
Expiration Date	12 months from date of receipt.

