

Product Datasheet

Vazyme - HiScript III All-in-one RT SuperMix Perfect for qPCR (R333-01)

Description	HiScript III All-in-one RT SuperMix Perfect for qPCR is an upgraded version of HiScript III RT Super...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Fast and convinient: The gDNA elimination and reverse transcription can be finished in one step. The whole reaction take only 15 min. Strong cDNA stability: The heat-labile DNase can be inactivated thoroughly, thus cDNA can be preserved for a longer period of time.
Expiration Date	12 months from date of receipt.