

Product Datasheet

Vazyme - ResiDNA Hunter Residual DNA Sample Preparation Kit [Out of stock] (RD101-01)

Description

ResiDNA Hunter The Residual DNA Sample Preparation Kit is applicable to the pretreatment of various ...

Storage

The Proteinase K, Glycogen and Yeast tRNA should be stored at -30 ~ -15°C and transported at ≤0°C. The other components should be stored at 15 ~ 25°C and transported at room temperature.

Note

For research use only

Expiration Date

Please enquire.

