

Product Datasheet

Canine adrenomedullin (ADM) ELISA Kit (orb1655750)

Description	Canine adrenomedullin (ADM) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.35-1000 pg/mL. Sensitivity: 4.42 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.