

Product Datasheet

Bovine Vasoactive intestinal peptide (VIP) ELISA Kit (orb1657611)

Description	Bovine Vasoactive intestinal peptide (VIP) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.17-500 pg/mL. Sensitivity: 2.63 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.