

Product Datasheet

General Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) ELISA Kit (orb1664610)

Description	General Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.055 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.