

Product Datasheet

Canine Estradiol,E2 ELISA Kit (orb1671152)

Description	Canine Estradiol,E2 ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 25 pg/mL-1000 pg/mL. Sensitivity: 25 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.