

Product Datasheet

Agarose Gel Extraction Kit - Column Kit (orb1733636)

Description

Agarose Gel Extraction Kit is designed for high-yield recovery of DNA from agarose gel with simultaneous removal of primer-dimers, primers, nucleotides, proteins, salt, agarose, ethidium bromide, and other impurities. The preparation is based on a silica-membrane technology for binding DNA in high-salt and elution in low-salt buffer. The kit provides a simple and efficient way to purify DNA in a size range between 100 bp and 10 kb. It requires no organic extractions or precipitation and guarantees high and stable recovery rates.

Storage

store at ambient temperature

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)