

Product Datasheet

SCRIPT Direct RT-qPCR ProbesMaster highROX (orb1733650)

Description	SCRIPT Direct RT-qPCR ProbesMaster highROX is designed for quantitative real-time analysis of target RNA directly from animal-or plant tissue, swabs and blood. The mix allows robust amplification avoiding the requirement of any prior RNA purification procedures. The mix is recommended for use with dual-labeled fluorescent probes, e.g. TaqMan, Molecular Beacon or Scorpion probes. It provides a powerful tool for multiplex-quantification of sample RNA in a broad dynamic range with exceptional sensitivity and precision. The mix contains all reagents required for RT-qPCR (except template, primer and labeled fluorescent probe) in a premixed 2x concentrated ready-to-use solution. High robustness, reliability and sensitivity of the mix are based on a genetically engineered reverse transcriptase and an antibody-blocked hot start polymerase in combination with an optimized and well-balanced buffer system. The mix ensures fast and easy preparation with a minimum of pipetting steps and is specially recommended for: Direct amplification of target RNA from various tissues samples; Direct detection of viral or bacterial RNA in nasal or throat swabs; Direct PCR from blood samples; Point-of-Care diagnostics.
Form/Appearance	liquid
Concentration	2x
Storage	store at -20 °C. avoid freeze/thaw cycles. stable at 4 °C for up to 4 weeks
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)