

Product Datasheet

Phosphatase Inhibitor Cocktail I (100× ddH₂O) (orb1909597)

Description	It is easy for the proteins to be degraded or dephosphorylated during in vitro extraction procedure, which could cause an inaccurate result of protein expression detection. Therefore, adding protease inhibitors or phosphatase inhibitors into the extracts would be an effective method to prevent degradation and dephosphorylation of protein. The phosphatase inhibitor mixture can effectively inhibit the dephosphorylation of common phosphatases on proteins and maintain the original phosphorylation state of proteins.
Storage	-20°C
Note	For research use only
Purity	98.00%
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)