

Product Datasheet

General Noradrenaline (NE) ELISA Kit (orb1952490)

Description	General Noradrenaline (NE) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.