

Product Datasheet

2X SDS-PAGE Sample Buffer without DTT or b-MeOH [Out of stock] (orb2064310)

Description	2X SDS-PAGE Sample Buffer without DTT or b-MeOH [Out of stock]
Conjugation	Unconjugated
Preservatives	Liquid (sterile filtered)
Form/Appearance	Liquid (sterile filtered)
Storage	Store container at -20°C or colder prior to opening. Aliquot and re-freeze for best results. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18 to 26°C) to restore solubility of some salts. Expiration date is six (6) months from date of opening container.
Note	For research use only
Expiration Date	12 months from date of receipt.