

Product Datasheet

Orange Gel Loading Buffer with DNA Stain protein (orb533362)

Description	Biorbyt Gel Loading Buffers with DNA Stain are formulated to facilitate loading of DNA samples into ...
Form/Appearance	Liquid
Concentration	6x conc
Storage	store at -20 °C
Note	For research use only
Application notes	Orange Gel Loading buffer with DNA Stain is recommended for analyzing DNA fragments ≤ 500 bp. Add one part of Gel Loading Buffer to 5 parts of the DNA sample as shown in the table below. Spin down the tubes and vortex gently to achieve homogeneity. DNA sample: 5 μ l 10 μ l 20 μ l 50 μ l; 6x Gel Loading Buffer: 1 μ l 2 μ l 4 μ l 10 μ l.
Expiration Date	6 months from date of receipt.