

## Product Datasheet

# Bacteria DNA Preparation - Column Kit (orb533448)

<b>Description</b>	The spin column based Bacteria DNA Preparation Kit is designed for rapid and high purity isolation o...
<b>Storage</b>	store at ambient temperature
<b>Note</b>	For research use only
<b>Application notes</b>	The obtained DNA is suitable for a variety of applications, including real-time PCR, southern blot analysis, genotyping and discovery or validation of SNP/SSR markers. Before start, add the following components (not included in the kit) as indicated on the respective bottle/tube: double distilled water to RNase A, Lysozyme and Proteinase K. 96-99% Ethanol to the Washing Buffer.
<b>Expiration Date</b>	12 months from date of receipt.