

AK5 rabbit pAb**Cat#: orb764494 (Manual)**


For research use only. Not intended for diagnostic use.

Product Name	AK5 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human AK5. AA range:391-440
Specificity	AK5 Polyclonal Antibody detects endogenous levels of AK5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Adenylate kinase isoenzyme 5
Gene Name	AK5
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal


Concentration	1 mg/ml
Observed band	25kD
Human Gene ID	26289
Human Swiss-Prot Number	Q9Y6K8
Alternative Names	AK5; Adenylate kinase isoenzyme 5; AK 5; ATP-AMP transphosphorylase 5

Background

adenylate kinase 5(AK5) Homo sapiens This gene encodes a member of the adenylate kinase family, which is involved in regulating the adenine nucleotide composition within a cell by catalyzing the reversible transfer of phosphate groups among adenine nucleotides. This member is related to the UMP/CMP kinase of several species. It is located in the cytosol and expressed exclusively in brain. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],



Western Blot analysis of various cells using AK5 Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using AK5 Antibody. The picture on the right is blocked with the synthesized peptide.



Explore. Bioreagents.

www.biorbyt.com