

TNNI3K rabbit pAb**Cat#: orb769371 (Manual)**


For research use only. Not intended for diagnostic use.

Product Name	TNNI3K rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human TNNI3K. AA range:301-350
Specificity	TNNI3K Polyclonal Antibody detects endogenous levels of TNNI3K protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide..
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Serine/threonine-protein kinase TNNI3K
Gene Name	TNNI3K
Cellular localization	Nucleus . Cytoplasm . Expressed at lower levels in the cytoplasm.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal

Concentration	1 mg/ml
Observed band	
Human Gene ID	51086
Human Swiss-Prot Number	Q59H18
Alternative Names	TNNI3K; CARK; Serine/threonine-protein kinase TNNI3K; Cardiac ankyrin repeat kinase; Cardiac troponin I-interacting kinase; TNNI3-interacting kinase
Background	This gene encodes a protein that belongs to the MAP kinase kinase kinase (MAPKKK) family of protein kinases. The protein contains ankyrin repeat, protein kinase and serine-rich domains and is thought to play a role in cardiac physiology. [provided by RefSeq, Sep 2012],



Western blot analysis of Mouse-kidney lysis using TNNI3K antibody. Antibody was diluted at 1:500. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TNNI3K Antibody. The picture on the right is blocked with the synthesized peptide.