

## Product Datasheet

# Vazyme - ClonExpress II One Step Cloning Kit (C112-02)

<b>Description</b>	The ClonExpress technology is a simple, fast, and highly efficient DNA seamless cloning technology. ...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C. ▲Please avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	Good performance: No ligase, no self-ligation, positive rate > 95% Broad Compatibility: Efficient cloning of 50 bp - 10 kb DNA fragments; applicable for single fragment Additional features: No need to consider the restriction site of insert fragments Software support: Provide online primer design software "CE Design"
<b>Expiration Date</b>	12 months from date of receipt.

