

## Product Datasheet

# Vazyme - ClonExpress MultiS One Step Cloning Kit (C113-02)

**Description** The ClonExpress technology is simple, fast, and highly efficient DNA seamless cloning technology. It...

**Storage** Store at -30 ~ -15°C and transport at ≤0°C. ▲Please avoid repeated freezing and thawing.

**Note** For research use only

**Application notes** Good performance: No ligase, no self-ligation, positive rate > 95% Broad Compatibility: Efficient cloning of 50 bp - 10 kb DNA fragments; applicable for 2 - 5 fragments Additional features: No need to consider the restriction site of insert fragments Software support: Provide online primer design software "CE Design"

**Expiration Date** 12 months from date of receipt.

