

Product Datasheet

Vazyme - ClonExpress Ultra One Step Cloning Kit V2 (C116-01)

Description The ClonExpress technology is simple, fast, and highly efficient DNA seamless cloning technology. It...

Storage Store at -30 ~ -15°C and transport at ≤0°C.

Note For research use only

Application notes Simple: Ready-to-use super mix in one tube Fast: Recombination only takes 5 - 30 minutes Efficient: Efficiently clone 50 bp - 10 kb fragments, the positive rate is > 95% Compatible: Suitable for GC-rich/AT-rich homologous sequences Special optimization: Compatible with 15 - 35 bp homologous sequences

Expiration Date 12 months from date of receipt.

