

## Product Datasheet

# Vazyme - ClonExpress Ultra One Step Cloning Kit V2 (C116-02)

**Description** The ClonExpress technology is simple, fast, and highly efficient DNA seamless cloning technology. It...

**Storage** Store at -30 ~ -15°C and transport at ≤0°C.

**Note** For research use only

**Application notes** Simple: Ready-to-use super mix in one tube Fast: Recombination only takes 5 - 30 minutes Efficient: Efficiently clone 50 bp - 10 kb fragments, the positive rate is > 95% Compatible: Suitable for GC-rich/AT-rich homologous sequences Special optimization: Compatible with 15 - 35 bp homologous sequences

**Expiration Date** 12 months from date of receipt.

