

## Product Datasheet

# Vazyme - Mut Express II Fast Mutagenesis Kit V2 (C214-02)

**Description** Mut Express II Fast Mutagenesis Kit V2 is a site-directed mutagenesis system based on the ClonExpres...

**Storage** Store at -30 ~ -15°C and transport at ≤0°C. ▲Please avoid repeated freezing and thawing.

**Note** For research use only

**Application notes** Low template usage: Robust and reliable exponential amplification method High-fidelity enzyme: Phanta Max Super-Fidelity DNA Polymerase (~128X Taq) minimizes unwanted secondary mutations Simple operation: Amplified product can be directly used in the recombination reaction after being digested with DpnI Widely applicable: Suitable for amplification of any site in the plasmid within 20 kb Software support: Provide online primer design software "CE Design"

**Expiration Date** Please enquire.

