

Product Datasheet

Vazyme - Mut Express MultiS Fast Mutagenesis Kit V2 (C215-01)

Description Mut Express MultiS Fast Mutagenesis Kit V2 is a kit based on ClonExpress rapid cloning technology, w...

Storage Store at -30 ~ -15°C and transport at ≤0°C. ▲Please avoid repeated freezing and thawing.

Note For research use only

Application notes Low template usage: Robust and reliable exponential amplification method High-fidelity enzyme: Phanta Max Super-Fidelity DNA Polymerase(~128X Taq) minimizes unwanted secondary mutations Simple operation: Amplified product can be directly used in the recombination reaction after being digested with DpnI Widely applicable: suitable for amplification of any site in the plasmid within 20 kb Software support: Provide online primer design software "CE Design"

Expiration Date 12 months from date of receipt.

