

Product Datasheet

Vazyme - 5min Universal Ligation Mix (C311-01)

Description	5min Universal Ligation Mix is a ready-to-use 2 x premix solution, containing T4 DNA ligase that ca...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Compatible: Compatible with sticky ends, blunt ends, TA , linker or adapter ligation Fast: Ligation can be completed quickly in 5 minutes at 25°C Convenient: Reaction mixture can be used directly for bacterial transformation
Expiration Date	12 months from date of receipt.

