

## Product Datasheet

# Vazyme - MycoBlue Mycoplasma Detector (D101-01)

<b>Description</b>	MycoBlue Mycoplasma Detector is designed for rapid detection of mycoplasma contamination in cell cul...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Fast: Incubate at 60°C for only 1 h to observe the results Simple: Direct detection without pretreatment Compatible: Cover most cell cultures common mycoplasma Precise: Accuracy comparable to qPCR, better than PCR
<b>Expiration Date</b>	12 months from date of receipt.

