

Product Datasheet

Vazyme - FastPure Host Removal and Microbiome DNA Isolation Kit [Out of stock] (DC501-01)

Description

This kit is applicable for selective removal of host nucleic acid from biological fluid samples (blo...

Storage

Store BOX 1 at -30 ~ -15°C and transport at ≤0°C. Store BOX 2 at 15 ~ 25°C and transport at room temperature. Store BOX 3 at 2 ~ 8°C and adjust the transportation method according to the destination

**Note**

For research use only

Expiration Date

Please enquire.