

## Product Datasheet

### Vazyme - RNase A (100 mg/ml) (DE111-01)

|                        |  |
|------------------------|--|
| <b>Description</b>     | RNase A is a high-quality endoribonuclease extracted from bovine pancreas that can specifically degr...            |
| <b>Storage</b>         | Store at 15 ~ 25°C for three months, and store at -30 ~ -15°C for long-term storage. Transport at room temperature |
| <b>Note</b>            | For research use only  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |