

Product Datasheet

Radox - Chloramphenicol (HRP9263)

| | |
|---------------------------|--|
| Description | Chloramphenicol - HRP Conjugate |
| Conjugation | HRP |
| Target | Chloramphenicol |
| Form/Appearance | Liquid |
| Storage | This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C. |
| Note | For research use only |
| Application notes | competitive ELISA with polyclonal antibody PAS9681/PAS9511 & PAS10031/PAS10032. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution. |
| Hazard Information | Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap. |
| Dilution Range | 1/640 & 1/1280. |
| Expiration Date | 12 months from date of receipt. |

