

Product Datasheet

Radox - Thiamphenicol/Florenicol (HRP9276)

Description	Thiamphenicol/Florenicol - HRP Conjugate
Conjugation	HRP
Target	Thiamphenicol/Florenicol
Form/Appearance	Liquid
Storage	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C.
Note	For research use only
Application notes	competitive ELISA with polyclonal antibody PAS9753/PAS9602. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
Hazard Information	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
Dilution Range	1/1280.
Expiration Date	12 months from date of receipt.

