

Product Datasheet

Radox - Streptomycin (HRP9338)

Description	Streptomycin - HRP Conjugate
Conjugation	HRP
Target	Streptomycin
Form/Appearance	Liquid
Storage	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C.
Note	For research use only
Application notes	The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
Hazard Information	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
Expiration Date	12 months from date of receipt.

