

Product Datasheet

Radox - Estriol (HRP9399)

Description	Estriol - HRP Conjugate
Conjugation	HRP
Target	Estriol
Form/Appearance	Liquid
Storage	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤ -20°C.
Note	For research use only
Application notes	competitive ELISA with polyclonal antibody PAS9880/PAS9881. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
Hazard Information	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
Dilution Range	1/1280.
Expiration Date	12 months from date of receipt.

