

Product Datasheet

Radox - Tacrolimus (HRP9425)

Description	Tacrolimus - HRP Conjugate
Conjugation	HRP
Target	Tacrolimus
Form/Appearance	Liquid
Storage	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C.
Note	For research use only
Application notes	competitive ELISA with polyclonal antibody PAS9875/PAS9874. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
Hazard Information	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
Dilution Range	1/320.
Expiration Date	12 months from date of receipt.

