

## Product Datasheet

### Radox - Chlorpromazine (HRP9457)

<b>Description</b>	Chlorpromazine - HRP Conjugate
<b>Conjugation</b>	HRP
<b>Target</b>	Chlorpromazine
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This conjugate may be stored for up to 3 months at +2 - +8°C. For long term storage, aliquot and store at ≤-20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	competitive ELISA with polyclonal antibody PAS10017/PAS10018. The conjugate is supplied as a concentrate. Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.
<b>Hazard Information</b>	Product should be protected from light exposure. Avoid repeated freeze/thaw cycles. Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.
<b>Dilution Range</b>	1/1520.
<b>Expiration Date</b>	12 months from date of receipt.

