

## Product Datasheet

### Radox - Apolipoprotein E4 (MAB10229)

<b>Description</b>	Apolipoprotein E4 Antibody (mAb)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Apolipoprotein E2 ( 1% ); Apolipoprotein E3 ( 1% ). Apolipoprotein E4 ( ~3ng/ml ).
<b>Preservatives</b>	0.09% Sodium Azide. Phosphate Buffered Saline (PBS)
<b>Form/Appearance</b>	Clear Liquid
<b>Concentration</b>	Lot dependent. Determined @ 280nm
<b>Storage</b>	Store at +2°C - +8°C. Refer to product label for lot-assigned expiration details. Product should be protected from light exposure.
<b>Note</b>	For research use only
<b>Application notes</b>	Product is suitable for ELISA and antibody recognition assays
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95%
<b>Expiration Date</b>	12 months from date of receipt.

