

Product Datasheet

Radox - Salbutamol (MAB9901)

Description	Salbutamol Antibody (mAb)
Species/Host	Sheep
Conjugation	Unconjugated
Immunogen	Salbutamol-BTG
Preservatives	0.10% Sodium Azide.. Phosphate Buffered Saline (PBS)
Form/Appearance	Clear Liquid
Concentration	Lot dependent. Determined @ 280nm with 1mg/ml = 1.4AU
Storage	Store at +2°C - +8°C. Refer to product label for lot-assigned expiration details. Product should be protected from light exposure.
Note	For research use only
Application notes	Product is suitable for ELISA and antibody recognition assays
Isotype	IgG
Clonality	Monoclonal
Purity	≥ 90%
Expiration Date	12 months from date of receipt.

