

Product Datasheet

Vazyme - 100 bp DNA Ladder (MD104-01)

Description	This product consists of double-stranded DNA fragments of specific molecular weight and is mixed wit...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Good stability: No change in product performance at room temperature for 6 months Reliable results: Electrophoretic bands are Sharp and clear Convenient: Provided with Gel Loading Buffer for tracking of sample DNA migration Precise: Exact molecular weight of DNA in each band
Expiration Date	12 months from date of receipt.

