

Product Datasheet

Vazyme - miRNA Universal SYBR qPCR Master Mix (MQ101-01)

Description miRNA Universal SYBR qPCR Master Mix is a special premix for qPCR assays based on SYBR Green I chime...

Storage Store at -30 ~ -15°C and protect from light. Transport at ≤0°C.

Note For research use only

Application notes Specificity: Chemically modified Hot-start AceTaq DNA Polymerase, completely blocks the enzyme activity below 95°C, and is equipped with the specificity-enhancing factor, which has stronger amplification specificity. Superior reaction system: The optimized Buffer composition and concentration are more suitable for the qPCR amplification of microRNA. Broad qPCR instrument compatibility: The special ROX in the mix can be used on almost all qPCR instruments without adjusting the concentration. Software support: Provide online primer design software 'miRNA Primer Design'.

Expiration Date 12 months from date of receipt.

