

## Product Datasheet

# Vazyme - VAHTS mRNA Capture Beads 2.0 (N403-01)

**Description** VAHTS mRNA Capture Beads 2.0 are 1  $\mu$ m Oligo(dT)-coupled paramagnetic microspheres intended for the ...

**Storage** Store at 2 ~ 8°C. Adjust the shipping method according to the destination.

**Note** For research use only

**Application notes** High purity: Specific enrichment of mRNA, higher purity and reduced rRNA residue. Wide compatibility: Compatible with multiple species - eukaryotes such as animals, plants, fungi. Compatible with a wider range of input samples—10 ng-12.5  $\mu$ g.

**Expiration Date** 12 months from date of receipt.

