

Product Datasheet

Vazyme - VAHTS AmpSeq Library Prep Kit V3 (NA210-01)

Description	VAHTS AmpSeq Library Prep Kit V3 is based on ultra-multiplex PCR, introducing several core technolog...
Storage	Store at -30 ~-15°C and transport at ≤ 0°C.
Note	For research use only
Application notes	High Data Uniformity: Uniformity > 95%. End Primer Digestion: Primer digestion reduces waste of sequencing capacity. Dual Platform Compatibility: Compatible with mainstream Illumina and Ion Torrent platforms. Broad Sample Type Compatibility: Compatible with gDNA, FFPE DNA, and cfDNA samples. Low Template Input: Starting template input as low as 1 ng.
Expiration Date	12 months from date of receipt.